


# IHE Integration Statement

July 2011



**Copyright © 2003-2011 medavis GmbH. All rights reserved.**

For feedback or questions please contact:

medavis GmbH  
Bannwaldallee 60  
D-76185 Karlsruhe  
Tel.: +49 721 92910 0  
Fax : +49 721 92910 99  
<http://www.medavis.com>  
[info@medavis.com](mailto:info@medavis.com)

This publication is protected by copyright. This document, as well as the software described in it, may be used only in accordance with the terms of the contract/license agreement.

The content of this document is intended for informational use only and is subject to change without notice. medavis Medizinische Informations Systeme GmbH assumes no responsibility or liability for any errors or inaccuracies that may appear in this document.

No part of this publication may be reproduced, stored in a retrieval system or transmitted, in any form or by any means, electronic, mechanical, recording, or otherwise, without the prior written permission of medavis Medizinische Informations Systeme GmbH.

medavis® is a registered trademark of medavis GmbH, Germany.

All trademarks and product names mentioned herein are trademarks or registered trademarks of their respective owners.

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
medavis GmbH	medavis RIS	3.0.5.1	May 24th, 2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles</b>	<b>Actors</b>	<b>Options</b>	
Access to Radiology Information	Image Manager/Archive	none	
Basic Image Review	Image Manager/Archive	none	
Basic Security	Secure Node (all actors)	none	
	Time Server	none	
	Audit Record Repository	none	
Cardiac Catheterization Workflow	ADT	none	
	Order Placer	none	
Charge Posting	Charge Processor	none	
	ADT	none	
	Order Filler	none	
Consistent Presentation of Images	Image Manager/Archive	none	
Echocardiography Workflow	ADT	none	
	Order Placer	none	
Evidence Documents	Image Manager/Archive	none	
Import Reconciliation Workflow	Order Filler	none	
	Importer	none	
	Patient Demographics Supplier	none	
Laboratory Scheduled Workflow	ADT	none	
	Order Placer	none	
	Order Result Tracker	none	
Laboratory Information Reconciliation	Order Placer	none	
	Order Result Tracker	none	
Mammography Image	Image Manager/Archive	none	
Nuclear Medicine Image	Image Manager/Archive	none	
Patient Administration Management	Patient Demographics Consumer	none	
	Patient Demographics Supplier	none	
	Patient Encounter Consumer	none	
	Patient Encounter Source	none	
Patient Information Reconciliation	ADT	none	
	Order Placer	none	
	Order Filler	none	
	PPS Manager	none	
	Image Manager/Archive	none	
Portable Data for Imaging	Portable Media Importer	none	
	Image Display	none	
	Display	none	
Patient Demographics Query	Patient Demographics Consumer	none	
	Patient Demographics Supplier	none	

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
medavis GmbH	medavis RIS	3.0.5.1	May 24th, 2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Radiation Exposure Monitoring	Image Manager/Archive	none	
Retrieve Information for Display	Information Source	none	
Scheduled Workflow	ADT	none	
	Order Placer	none	
	Order Filler	Billing & Material Management Appointment Notification	
	PPS Manager	PPS Exception Management	
	Image Manager/Archive	none	
	Evidence Creator	none	
<b>Internet address for vendor's IHE information:</b>			
<a href="http://www.medavis.com/index.php/products/module/ihe-interface.html">www.medavis.com/index.php/products/module/ihe-interface.html</a>			
Links to Standards Conformance Statements for the Implementation			
<b>HL7</b>	<a href="http://www.medavis.com/down/CS_medavis_HL7_Interface_Specification_V1.7_gb.pdf">www.medavis.com/down/CS_medavis_HL7_Interface_Specification_V1.7_gb.pdf</a>		
<b>DICOM</b>	<a href="http://www.medavis.com/index.php/downloads.html">www.medavis.com/index.php/downloads.html</a>		
<b>Links to general information on IHE</b>			
In North America:	<a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe:	<a href="http://www.ihe-europe.org">www.ihe-europe.org</a>
		In Japan:	<a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
medavis GmbH	medavis PACS	4.4	May 24th, 2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles</b>	<b>Actors</b>	<b>Options</b>	
Access to Radiology Information	Image Display	none	
	Image Manager/Archive	none	
Audit Trail and Node Authentication	Secure Application	none	
Cardiac Catheterization Workflow	Image Display	none	
	Image Manager/Archive	none	
	PPS Manager	none	
Consistent Presentation of Images	Image Display	none	
	Image Manager/Archive	none	
	Print Composer	none	
Cross-Enterprise Document Sharing.b (XDS.b)	Document Repository	none	
	Document Source	none	
	Embedded Repository	none	
Cross-Enterprise Document Sharing for Imaging.b (XDS-I.b)	Document Repository	none	
	Imaging Document Source	none	
Cross-Enterprise Document Reliable Interchange	Document Recipient	none	
	Document Source	none	
Echocardiography Workflow	Image Manager/Archive	none	
	PPS Manager	none	
Evidence Documents	Image Display	none	
	Image Manager/Archive	none	
Eye Care Evidence Documents	Image Manager/Archive	none	
Eye Care Workflow	Image Manager/Archive	none	
Image Fusion	Image Manager/Archive	none	
Import Reconciliation Workflow	Image Manager/Archive	none	
	Importer	none	
	Patient Demographics Supplier	none	
	PPS Manager	none	
Key Image Note	Image Display	none	
	Image Manager/Archive	none	
Mammography Acquisition Workflow	Image Display	none	
	Image Manager/Archive	none	
Mammography Image	Image Display	none	
	Image Manager/Archive	none	
Nuclear Medicine Image	Image Display	none	
	Image Manager/Archive	none	
Patient Demographics Query	Patient Demographics Consumer	none	
Patient Information Reconciliation	Image Manager/Archive	none	
	PPS Manager	none	

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
medavis GmbH	medavis PACS	4.4	May 24th, 2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Portable Data for Imaging	Image Display	none	
	Portable Media Creator	none	
	Portable Media Importer	none	
Post Processing Workflow	Image Manager/Archive	none	
Presentation of Grouped Procedures	Image Manager/Archive	none	
	PPS Manager	none	
Radiation Exposure Monitoring	Image Manager/Archive	none	
Reporting Workflow	Image Manager/Archive	none	
Retrieve ECG for Display	Display	none	
	Information Source	none	
Scheduled Workflow	Image Display	none	
	Image Manager/Archive	none	
	PPS Manager	none	
<b>Internet address for vendor's IHE information:</b>			
<a href="http://www.medavis.com/index.php/products/module/ihe-interface.html">www.medavis.com/index.php/products/module/ihe-interface.html</a>			
Links to Standards Conformance Statements for the Implementation			
<b>HL7</b>	<a href="http://www.medavis.com/down/CS_medavis_HL7_Interface_Specification_V1.7_gb.pdf">www.medavis.com/down/CS_medavis_HL7_Interface_Specification_V1.7_gb.pdf</a>		
<b>DICOM</b>	<a href="http://www.medavis.com/index.php/downloads.html">www.medavis.com/index.php/downloads.html</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	

www.medavis.de + +  
 + + + + + +  
 www.medavis.ch + +  
 + + + + + +  
 www.medavis.at + +  
 + + + + + +  
 + + + + + +  
 + + + + + +  
 + + + + + +  
 + + + + + +  
 + + + + + +

## medavis at a glance

medavis GmbH, an owner-managed company, was founded in 1997, after completing the basic medavis RIS product.

We combine practical experience with new software technology and offer a complete range of services from development, installation and training right through to ongoing customer care.

Our independence and expertise are the basis for successful links to other healthcare system vendors of every size. Customers in Germany, Austria and Switzerland benefit as a result.

medavis GmbH  
 Bannwaldallee 60  
 76185 Karlsruhe  
 Germany

Tel +49 721 92910-0  
 info@medavis.com

medavis 

medavis and other medavis names referred to herein are trademarks or registered trademarks of medavis GmbH

All other product and company names mentioned herein are trademarks of their respective owners. All other rights not expressly referred to herein are reserved.  
 Printed on chlorine and acid-free paper.